



BRAND Acrylic and Polystyrene Spectrophotometry Cuvettes

BRAND plastic cuvettes are an excellent choice ([see why](#)) for most spectrophotometry applications in the visible range. Manufactured from polystyrene (PS) or polymethyl methacrylate (PMMA or 'acrylic'), they provide accurate, reliable results, even at wavelengths as low as 300nm. ([For wavelengths below 300nm, see our UV-Cuvettes.](#)) These disposable cuvettes are available in macro and semi-micro sizes to accommodate most sample volumes. BRAND plastic cuvettes are compatible with most spectrophotometers and photometers.



Easy to Use

BRAND plastic disposable cuvettes include an arrow mark to indicate direction of transmission and reduce variation.

Options for Most Needs

PS and PMMA cuvettes are available in both macro and semi-micro sizes to adapt to the needs of most laboratory applications. Semi-micro cuvettes require only 1.5mL of sample volume to provide reliable, accurate results.

Performs Consistently

Cuvettes are manufactured from high-quality materials, and are grouped by manufacturing mold cavity number to ensure the lowest variation in extinction coefficient.

Manufactured by BRAND

One of the world's largest producers of cuvettes, BRAND's 15 years of experience manufacturing quality cuvettes is your assurance of reliability.

Ordering Information:

BRAND PS and PMMA cuvettes are available from leading laboratory dealers. For more information, or to request a quote, contact your favorite laboratory supply dealer or call BrandTech Scientific at 888-522-2726.

Description	Catalog Number
BRAND Polystyrene (PS) Cuvettes	
Cuvette, PS, Macro, pack of 500	759070D
Cuvette, PS, Macro, pack of 100	759071D
Cuvette, PS, Semi-Micro, pack of 500	759075D
Cuvette, PS, Semi-Micro, pack of 100	759076D
BRAND Methacrylate(PMMA) Cuvettes	
Cuvette, PMMA, Macro, pack of 500	759080D
Cuvette, PMMA, Macro, pack of 100	759081D
Cuvette, PMMA, Semi-Micro, pack of 500	759085D
Cuvette, PMMA, Semi-Micro, pack of 100	759086D
Accessories	
Cuvette Rack, PP, w/ 16 numbered positions	759500