



## accu-jet<sup>®</sup> *pro* Pipette Controller

The accu-jet<sup>®</sup> *pro* elevates pipette control to the next generation. It improves upon the superior function of the accu-jet with new features, and improved functionality.

**PROgressive Speed Control!** A new range dial lets you set maximum pipetting speed, but still provides the touch-sensitive, continuously adjustable speed control that made the accu-jet so popular.

**PROportioned Perfectly!** The accu-jet *pro* has a reshaped ergonomic grip that works well with both small and large hands.

**PROlonged Battery Life!** The accu-jet *pro* provides 8 hours of full speed operation on a 4-hour charge. There's even an LED indicator when you have two hours of pipetting time left. Convenient battery access with a slide-on cover.

**PROtects Cell Pellets.** Unparalleled soft blow-out control minimizes disturbance of cell pellets, or density gradient interfaces.

All this, along with the preferred features of the original accu-jet; like a true 0.2µm hydrophobic filter for culture work, a check valve to help protect the internal components, exceptional balance and light weight, and external venting of pipetted vapors.



### Ordering Information:

The accu-jet<sup>®</sup> *pro* comes complete with NiHM battery pack, wall holder, two sterile replacement 0.2µm hydrophobic filters, 120V/60Hz charger, operating manual, and one-year warranty.

For more information, or to request a quote, contact your favorite laboratory supply dealer or call BrandTech Scientific at 888-522-2726.

Description	Catalog Number
<b>accu-jet<sup>®</sup> <i>pro</i> Pipette Controller</b>	
accu-jet pro Pipette Controller, Dark Blue, 120V	2026330
accu-jet pro Pipette Controller, Magenta, 120V	2026331
accu-jet pro Pipette Controller, Green, 120V	2026332
accu-jet pro Pipette Controller, Royal Blue, 120V	2026333
<b>Accessories</b>	
Membrane filter, 0.2um, Sterile, pack of 10	2026539
Membrane filter for accu-jet pro, Non-Sterile, 0.2um, pack of 10	26535
Membrane filter, 3um, pack of 10	26056
Membrane filter, 3um, pack of 1	26052
Battery charger for accu-jet pro, 120V, US plug, pack of 1	26603

Adapter, silicone, with non-return valve, for accu-jet pro, pack of 1	26508
Wall Holder for accu-jet pro, pack of 1	26540
Snap Tapes, 12 x 65mm, 1 pair, for accu-jet pro wall holder, pack of 1	26542
NiMH battery for accu-jet pro, pack of 1	26630